



Quality
You Can
Trust
Style
You'll Love

FOUR
SQUARE
ADVT
CREATIVITY BEYOND IMAGINATION

MEDICAL FURNITURE AND HOSPITAL EQUIPMENT

ABOUT US

At SHREYA MEDICAL SYSTEM, we specialize in providing high-quality medical furniture designed to enhance the functionality and comfort of healthcare environments. With over 12 years of experience in the industry, we understand the unique needs of medical professionals and patients alike.

Our mission is to support the healthcare community by offering innovative solutions that improve patient care and streamline operations. From examination tables to waiting room seating, each piece of furniture is crafted with durability, safety, and ease of use in mind.

We are committed to using sustainable materials and advanced manufacturing processes to ensure that our products not only meet the highest standards of quality but also contribute to a healthier planet. Our dedicated team works closely with healthcare providers to design furniture that meets the specific requirements of various medical settings.

Join us in creating spaces where healing and care come together seamlessly. Explore our range of medical furniture and discover how we can help elevate your healthcare environment.

COMPANY SHOWROOM



QUALITY CONTROL AND R&D

we understand that in the medical field, quality is paramount. Our commitment to excellence is reflected in our rigorous Quality Control (QC) processes and our dedication to Research and Development (R&D).

FAIR BUSINESS PRACTICES

Transparent Pricing: We believe in clear, upfront pricing without hidden fees. Our clients can trust that they are receiving the best value for their investment.

Quality Assurance: We prioritize quality to ensure safety and satisfaction for healthcare providers and their patients.

Continuous Improvement: We are dedicated to ongoing evaluation and enhancement of our practices to better serve our clients and contribute positively to the healthcare industry.

ENVIRONMENT

By choosing eco-friendly materials and minimizing waste, we are working toward a healthier planet as we contribute to healthier lives.

VISION

Our vision is to elevate the patient experience by creating furniture that prioritizes comfort, dignity, and ease of use.

VALUES

Our furniture is built with high-quality materials and reliability. Innovation is at the heart of our approach. We believe that patient experience is crucial to effective care. Our furniture promotes patient comfort.





TURNKEY SOLUTIONS

STEP
01

**MEDICAL
FURNITURE**

STEP
02

**MODULAR
OPERATION ROOMS**

STEP
03

**MEDICAL GAS
PIPE LINE**

STEP
04

**SURGICAL
TABLES & LIGHTS**

STEP
05


**NEONATAL
CARE PRODUCTS**



**FOUR
SQUARE
ADVT**
CREATIVITY BEYOND IMAGINATION...



ICU BEDS



ICU BED MOTORISED (7 Function)

SMS-101

Description

- Positions Operated by Nurse Control Panel (NCP) & Patient Hand Remote: Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg, Hi – Low, Cardiac Chair Position, CPR.
- Dual Side Manually Operated CPR (For Backrest).
- One Touch button for Chair Position, Auto CPR (For Zero Position).
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails With Inbuilt Bed Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on ABS cover base with 125mm Central Locking Twin Wheel Castors along with directional Lock.
- Heavy Duty I.V Rod and Provision at 4

Locations.


- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2230L x 1010W x 590 – 870H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- X-ray Translucent Top for Backrest with X-ray Cassette Holder Beneath Backrest
- Polymer 4 Section Platform Tops.
- Extension frame for Bed.
- Battery Backup
- LED Night Lamp and Battery Indicator.




X-Ray Preamble Top



Nurse Control System



ICU BEDS



ICU BED MOTORISED (7 Function)

SMS-102

Description

- Positions Operated by Nurse Control Panel (inbuilt on bed panel) & Patient Hand Remote; Backrest, Kneerest, Trendelenburg & Reverse Trendelenburg, Hi-Low, Cardiac chair position, CPR.
- Inbuilt Control System in Side Railing
- Dual Side Manually Operated CPR (For Backrest).
- One Touch button for Chair Position, Auto CPR (For Zero Position).
- Battery Backup
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails With Inbuilt Bed Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on ABS cover base 125mm Central Locking Twin Wheel Castors along with directional Lock.
- Linen Holder for fresh Bedsheet under the

footside panel.

- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2230L x 1050W x 450 – 750H mm
- Angle of Back Section: 70°
- Angle of Knee Section: 25°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- X-ray Translucent Top for Backrest with X-ray Cassette Holder Beneath Backrest
- ABS 4 Section Platform Tops.
- Bed Extension.
- LED Night Lamp and Battery Indicator.



Patient Control System



Trolley ABS Base Cover



Nurse Control System



Handset

ICU BEDS



ICU BED MOTORISED (5 Function)

SMS-103

Description

- Backrest, Kneerest, Trendelenburg, Reverse Trendelenburg & Hi-Low Operated By Hand Remote System.
- Dual Side Manually Operated CPR (For Backrest).
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails With Inbuilt Bed Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on ABS cover base 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.


- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2230L x 1050W x 450 – 750H mm
- Angle of Back Section: 75°
- Angle of Knee Section: 25°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock.
- Battery Backup
- ABS Platform on top for bed.



ICU BED SEMI MOTORISED (3 Motorised + 2 Manual Functions)

SMS-104

Description

- Backrest, Kneerest, Hi-Low & CPR Operated By Hand Remote System.
- Trendelenburg, Reverse Trendelenburg operated by manual crank system.
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails With Inbuilt Bed Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2230L x 1080W x 580 – 880H mm
- Angle of Back Section: 75°
- Angle of Knee Section: 25°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock.
- Battery Backup
- Polymer Platform on top for Bed

ICU BEDS



ICU BED SEMI MOTORISED (1 Motorised + 3 Functions Manual)

SMS-105

Description

- Backrest Operated by Hand Remote
- Kneerest, Trendelenburg & Reverse Trendelenburg & Hi-Low operated by manual crank System.
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails with Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.

- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2270L x 1080W x 580 – 880H mm
- Angle of Back Section: 75°
- Angle of Knee Section: 35°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock



ICU BED SEMI-MOTORISED (2 Motorised + 3 Manual Functions)

SMS-106

Description

- Two Function Manual other Function are Motorised.
- Dual Side Manually Operated CPR (For Backrest).
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails With Inbuilt Bed Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on ABS cover base 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.

- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2230L x 1050W x 450 – 750H mm
- Angle of Back Section: 75°
- Angle of Knee Section: 25°
- Angle of Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock.
- Battery Backup
- ABS Platform on top for bed.

ICU BEDS



ICU BED MECHANICAL (5 FUNCTION)

SMS-107

Description

- Backrest Kneerest, Trendelenburg & Reverse Trendelenburg & Hi-Low operated by manual crank System.
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Four individual Side Rails with Angle Indicator
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 125mm Diagonal locking Wheel Twin Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2270L x 1080W x 580 – 880H mm
- Back Section: 0-75°
- Knee Section: 0-35°
- Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock.

ICU BEDS



ICU BED MECHANICAL (ABS Panel & Collapsible Railing)

SMS-108

Description

- Backrest Kneerest, Trendelenburg & Reverse Trendelenburg & Hi-Low operated by manual crank System.
- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Collapsible Safety Side Railing.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2120L x 980W x 580 – 880H mm
- Back Section: 0-75°
- Knee Section: 0-35°
- Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock.

ICU BEDS



ICU BED MECHANICAL (DELUXE)

SMS-109

Description

- Backrest Kneerest, Trendelenburg & Reverse Trendelenburg & Hi-Low operated by manual crank System.
- Detachable heavy ABS Head and Foot Panels easily removable without using any tools with corner buffer.
- Collapsible Safety Side Railing.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2120L x 980W x 580 – 880H mm
- Back Section: 0-75°
- Knee Section: 0-35°
- Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock

ICU BEDS



ICU BED MECHANICAL (Standard SS Panel)

SMS-110

Description

- Backrest Kneerest, Trendelenburg & Reverse Trendelenburg & Hi-Low operated by manual crank System.
- Detachable SS Head and Foot Panels with easily removable without using any tools with corner buffer.
- Collapsible Safety Side Railing.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 125mm Diagonal locking Twin Wheel Castors
- Heavy Duty I.V Rod and Provision at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2120L x 980W x 580 – 880H mm
- Back Section: 0-75°
- Knee Section: 0-35°
- Trendelenburg / Reverse Trendelenburg: 12° / 12°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Central Locking Castors with direction Lock

FOWLER BEDS



FOWLER BED MOTORISED (ABS Panel & Railing)

SMS-111

Description

- Backrest, Kneerest, operated by Remote Control System.
- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Polymer Moulded Four Pieces individual Side Rails With Inbuilt Angle Indicator.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on trolley type base with 125mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2200L x 1080W x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.

Also available Collapsible Railing

FOWLER BEDS



FOWLER BED MECHANICAL (2 Function.) ABS Panel & Railing

SMS-112

Description

- Backrest, Kneerest, operated by Manual Crank System .
- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Polymer Moulded Four Pieces individual Side Rails With Inbuilt Angle Indicator.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on trolley type base with 125mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2270L x 1080W x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.

FOWLER BEDS



FOWLER BED ELECTRICAL (ABS Panel & Collapsible Railing)

SMS-113

Description

- Backrest, Kneerest, operated by Remote Control System.
- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Collapsible Safety Side Railing
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted 125mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2200L x 1080W x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover.
- Heavy Duty I.V Rod.

Also available in Mechanical

FOWLER BEDS



FOWLER BED MECHANICAL (ABS Panel)

SMS-114

Description

- Backrest, Kneerest, Leg rest operated by Manual Crank System .
- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 100mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2130L x 980 x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section Foam Mattress with Zipper Rexine Cover.
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.

FOWLER BEDS



FOWLER BED MECHANICAL (SS Laminated Panel)

SMS-115

Description

- Backrest, Kneerest, Legrest operated by Manual Crank System.
- Detachable S.S with laminated Head and Foot Panels.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 100mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2070L x 920W x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.

FOWLER BEDS



FOWLER BED MECHANICAL (Std.)

SMS-116

Description

- Backrest, Kneerest, Legrest operated by Manual Crank System.
- Detachable Polymer Head and Foot Portable Panels.
- Framework made-up of CRCA Tubes.
- 4 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on 100mm Diagonal Locking Castors.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2070L x 920W x 600H mm.
- Angle of Back Section: 75°
- Angle of Knee Section: 35°

Optional Accessories

- 4 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.

SEMI-FOWLER BEDS



SEMI-FOWLER BED MOTORISED

SMS-117

Description

- Backrest, operated by Hand Remote System .
- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Framework made-up of CRCA Tubes.
- 2 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on PVC Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2110L x 980W x 600H mm.
- Angle of Back Section: 85°

Optional Accessories

- 2 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails
- 100mm Diagonal locking Castors

SEMI-FOWLER BEDS



SEMI-FOWLER BED MECHANICAL (ABS Panel)

SMS-118

Description

- Polymer Moulded Head and Foot Panels easily removable without using any tools
- Framework made-up of CRCA Tubes.
- 2 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on PVC Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Backrest operated by Smooth Crank System with Foldable Handle.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2130L x 980W x 600H mm.
- Angle of Back Section: 85°

Optional Accessories

- 2 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.
- 100mm Diagonal locking Castors.

SEMI-FOWLER BEDS



SEMI-FOWLER BED MECHANICAL (SS Laminated Panel)

SMS-119

Description

- Detachable S.S with laminated Head and Foot Panels easily removable.
- Framework made-up of CRCA Tubes.
- 2 Section Platform Top made-up of Perforated CRCA Sheet.
- Bed Mounted on PVC Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Backrest operated by Smooth Crank System with Foldable Handle.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2070L x 920W x 600H mm.
- Angle of Back Section: 85°

Optional Accessories

- 2 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.
- 100mm Diagonal locking Castors.

SEMI-FOWLER BEDS



SEMI-FOWLER BED MECHANICAL (STD.)

SMS-120

Description

- Polymer Easily Detachable Head and Foot Portable Panels.
- Framework made-up of CRCA Tubes.
- 2 Section Platform Top made-up of Oval Perforated CRCA Sheet.
- Bed Mounted on PVC Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Backrest operated by Smooth Crank System with Foldable Handle.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 2070L x 920W x 600H mm.
- Angle of Back Section: 85°

Optional Accessories

- 2 Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Collapsible Safety Side-Rails.
- 100mm Diagonal locking Castors.

PLAIN BEDS



Plain Bed (Std.)

SMS-121

Description

- CRCA Fix Head and Foot Panels.
- Single Platform Top made-up of Perforated CRCA Sheet .
- Bed Top made-up of single platform Perforated CRCA Sheet.
- Framework made-up of CRCA Tubes.
- Bed Mounted on PU Heavy Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Urine Bag Holder on Both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Approx. Size: 2010L x 920W x 600H MM

Optional Accessories

- Single Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- Swing away Safety Side-Rails.
- 100mm Diagonal locking Castors

PEDIATRIC BED (STD.)

SMS-122

Description

- Bed Top made-up of single platform Perforated CRCA Sheet.
- Framework made-up of CRCA Tubes.
- Bed Mounted on PU Heavy Stumps/ Shoes.
- Provision for IV Rod at 4 Locations.
- Full Length Swing Down Safety Side Rails both Sides.
- Finish: Pre-treated & Epoxy Powder Coated.

Dimensions (Approx.)

- Overall Size: 1485L x 860W x 600H mm

Optional Accessories

- Single Section PU Foam Mattress with Zipper Rexine Cover
- Heavy Duty I.V Rod.
- 100mm Diagonal locking Castors



ATTENDANT BED

SMS-123

Description

- Framework made of Rectangular / Square M.S. Tube.
- 50 mm (2") Foam Padded Fixed Mattress Covered with Good Quality Rexine.
- Legs fitted with PU Heavy stumps / Shoes.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 1840L x 620W x 490H MM



ATTENDANT BED CUM CHAIR

SMS-124

Feature and Benefits

- Frame is made of CRCA tubular pipe mounted on PVC stumps.
- Padded with good quality mattress.
- Metal body is pretreated and epoxy powder coated.
- dimension is 1900 x 620 x 580 mm.




ATTENDANT BED WITH 2 STORAGE BOX

SMS-125

Feature and Benefits

- Frame is made of CRCA tubular pipe and sheets mounted on PVC stumps.
- Padded with good quality mattress.
- Metal body is pretreated and epoxy powder coated.
- Dimension is 1870 x 600 x 460 mm.



BED SIDE LOCKER

BED SIDE LOCKER (ABS Top & Cabinet)

SMS-201

Description

- Locker made of CRCA sheets/ Square CRCA Pipe.
- Top, Drawer & Cabinet door made-up of Premium quality of Polymer.
- Mounted on 50 mm wheels at rear and 2 pcs PVC stumps/ Shoes in front.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 400L x 400W x 850H mm



BED SIDE LOCKER (Fully ABS)

SMS-202

Description

- Locker made of premium quality of Polymer.
- Provide with single Cabinet, Drawer, Tray & Towel rail with Polymer.
- Mounted on 50 mm 4 wheels.
- Overall Size: 450L x 420W x 775H mm.

BED SIDE LOCKER (Supreme)

SMS-202-B

Description

- Locker made of CRCA sheets.
- Provided one Cabinet with one Drawer.
- Provided with Stainless Steel Top.
- Mounted on 50mm twin wheel castor.
- Finish: Pre-treated & Epoxy Powder Coated.
- Overall Size: 400 x 450 x 850 mm



BED SIDE LOCKER (DLX.)

SMS-203

Description

- Locker made of CRCA sheets/ Square CRCA Pipe.
- Top made-up of Stainless Steel & Drawer & Cabinet door CRCA Sheet.
- Mounted on 50 mm wheels at rear and 2 pcs PVC stumps/ Shoes in front.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 400L x 400W x 850H mm

BED SIDE LOCKER (Supreme Plus)

SMS-203-B

Description

- Locker made of CRCA sheets.
- Provided one Cabinet with one Drawer.
- Provided with Side Storage.
- Provided with Stainless Steel Top.
- Mounted on 50mm twin wheel castor.
- Finish: Pre-treated & Epoxy Powder Coated.
- Overall Size: 400 x 450 x 850 mm



BED SIDE LOCKER (STD.)

SMS-204

Description

- Made-up of CRCA Sheet & Squar Pipe with Stainless Steel Top
- Provided with Single Storage Box.
- Legs Mounted on PVC stumps / Rubber shoe.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 400L x 400W x 850H mm

BED SIDE LOCKER (Fully SS)

SMS-204-B

Description

- Locker made of SS sheets & SS Pipe.
- Provided with Three shelves.
- Mounted on PVC stumps
- Overall Size: 400 x 400 x 810 mm



OVER BED TABLE

OVER BED TABLE (PNEUMATIC)

SMS-205

Description

- Polymer premium quality Top with CRCA Heavy quality frame.
- CRCA Square Pipe telescopic section fitted with gas spring pump for height adjustment .
- M.S Heavy base for better stability.
- Mounted on 50 mm wheels.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 890L x 435W x 790 – 1070H mm.
- Top Size: 880L X 435W X 28H mm.



Also available in single top

OVER BED TABLE (GEAR)

SMS-206

Description

- Two Section Wooden premium quality Top with CRCA Heavy quality frame.
- CRCA Square Pipe telescopic section fitted with smooth gear mechanism for height adjustment .
- M.S Heavy base for better stability.
- Mounted on 50 mm wheels.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 1000L x 435W x 825 – 1180H mm
- Top Size: 880L X 435W X 28H mm

OVER BED TABLE (KNOB)

SMS-207

Description

- Wooden laminated Top with CRCA Heavy quality frame.
- CRCA Square Pipe telescopic section fitted with knob for height adjustment .
- M.S Heavy base for better stability.
- Mounted on 50 mm wheels.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 1000L x 435W x 825 – 1180H mm
- Top Size: 880L X 435W X 28H mm



OVER BED TABLE (FIXED)

SMS-208

Description

- Wooden laminated Top with CRCA Heavy frame.
- Framework made of CRCA tubular Mounted on 50 mm wheels.
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size: 1060L x 435W x 1060H mm.
- Top Size: 880L X 435W X 28H mm

WARD / OT EQUIPMENTS

ABS CRASH CART TROLLEY

SMS-601

Description

- Framework & Drawers made of ABS body with aluminum columns.
- Trolley having Utility Container, File Holder, Defibrillator Shelf, oxygen tank at rear, CPR board at rear & 3 side rail on top. IV Rod with 4 hooks, Waste bin (2pcs), Needle disposal holder.
- Trolley having Sliding Side Shelf for writing purpose.
- Electric point with 2 socket backside with hook
- Centralized locking system for 2 small & 3 big drawer
- Provided with adjustable inner partitions.
- 125 mm Diagonal locking castors.
- Overall Size (approx.): 750L x 475W x 950H mm

MODEL 2: AT-7501182, Size: 750*475*950mm



INSTRUMENT CABINET

SMS-602

Description

- Made of CRCA sheets provided with single door and glass shelves.
- Pre-treated and epoxy powder coated.
- Overall Size (approx.): 1500 L x 600 W x 400 D mm



TRANSPORT TROLLEY

EMERGENCY & RECOVERY TROLLEY (HYDRAULIC)

SMS-301

Description

- Frame made of rectangular CRCA Pipe with Stainless Steel swing away railing.
- Operated by Smooth Pneumatic pump for Trendelenburg and Reverse Trendelenburg positions.
- Height Adjustable by hydraulic pump system.
- Removable two section X-Ray permeable top with Head Raise on Ratchet Mechanism.
- Trolley Mounted on 125 mm diagonal locking castors
- Provided with Oxygen Cylinder Cage, Sliding X-ray Cassette holder, Drop down safety side railings, Telescopic I.V. Rod, Utility Tray and Suitable Mattress of 50 mm.
- Finish: Pretreated and Epoxy powder coated.

Dimensions (Approx.)

- Overall Size: 2100L x 570W x 650 – 950H mm.
- Stretcher top Size: 1940L x 605W mm
- Back Section: 70°



EMERGENCY & RECOVERY TROLLEY (MANUAL)

SMS-302

Description

- Frame made of rectangular CRCA Pipe with Stainless Steel swing away railing.
- Operated by Smooth manual crank system for Trendelenburg and Reverse Trendelenburg positions.
- Height Adjustable by manual crank system.
- Removable two section X-Ray permeable top with Head Raise on Ratchet Mechanism.
- Trolley Mounted on 125 mm diagonal locking castors
- Provided with Oxygen Cylinder Cage, Drop down safety side railings, Telescopic I.V. Rod, Utility Tray and Suitable Mattress of 50 mm.
- Finish: Pretreated and Epoxy powder coated.

Dimensions (Approx.)

- Overall Size: 2100L x 570W x 650 – 950H mm.
- Stretcher top Size: 1940L x 605W mm
- Back Section: 70°



STRETCHER TROLLEY 2 Section

SMS-303

Description

- Framework made of Tubular CRCA Pipe with S.S. handle.
- 2 Section Platform Top with Fixed 50mm Mattress.
- Trolley mount on 150 mm diagonal Locking castors.
- Swing Away Type Safety Side Railing in M.S.
- Provided with Adjustable S.S. I.V. Rod, Utility Tray, and Cylinder Cage.
- Finish: Pretreated and Epoxy powder coated.
- Stretcher Top Size (Handle to handle) : 2090L x 560W
- Overall Size: 2090L x 590W x 840H mm.



Note : Also available in SS

STRETCHER ON TROLLEY (STD.)

SMS-304

Description

- Removable Stretcher top on SS Tubular Steel trolley mounted on 150mm Diagonal Locking Castors.
- Framework made of Stainless steel pipe and Sheet.
- Trolley overall Size (Handle to handle) : 2100L x 560W x 810H mm.

Optional Accessories

- S.S. Telescopic Heavy Duty I.V. Rod.
- Safety Side-Rails, Swing Away .
- Mattress 50mm with good quality raxine cover.

Note : Also available in epoxy powder coated.



Powder Coated

TRANSPORT TROLLEY

EMERGENCY RECOVERY TROLLEY DELUXE (4 SECTION TOP)

SMS-305

Description

- Frame made of rectangular CRCA pipe with collapsible side railing
- Removable radiolucent top with X-Ray cassette holder
- Four section top with back rest and manual maneuver for Knee rest via crank
- Trendlenburg/Reverse Trendlenburg position maneuverable using crank mechanism
- High-Low maneuverable by hydraulic pump
- Castors of 125mm dia with diagonal locking are provided for mobility with control.
- SS IV stand.
- The complete metal components are pre treated and epoxy powder coated
- Frame size : 2050L x 760W x 650-900H mm



STRETCHER ON TROLLEY (ABS RAILING)

SMS-306

Description

- Frame is made of CRCA Sheet & Pipes.
- Trolley is mounted over 125 mm diagonally locking castor wheels.
- Provision for saline stand.
- Provided with swing down polymer moulded safety side rails.
- Backrest can be adjusted via gas spring up to 65°
- High-Low can be adjusted via crank mechanism
- Dimension : 2050x760 x 560-800 mm
- the complete metal component is pre treated and epoxy powder coated.



EXAMINATION TABLE



EXAMINATION COUCH (with Drawers & Cabinet)

SMS-401

Description

- Fixed 50mm Foam Mattress with cover top in two section
- Body framework made up of CRCA sheet & mounted on CRCA leg covered with PU stumps/shoes.
- Smooth Pneumatic mechanism to raise Headrest.
- Provided 3 sliding drawers & 3 separate cabinet with doors.
- B.P. apparatus tray on head side.
- Sliding footstep located at bottom center.
- Finish: Pre-treated & epoxy powder coated
- Overall Size (Approx.): 1870L x 560W x 820H mm

GYNAE EXAMINATION COUCH (with Drawers & Cabinet)

SMS-402

Description

- Fixed 50mm Foam Mattress with cover top in two section
- Body framework made up of CRCA sheet & mounted on CRCA leg covered with PU stumps/shoes.
- Smooth Pneumatic mechanism to raise Headrest.
- Provided 3 sliding drawers & 3 separate cabinet with doors.
- B.P. apparatus tray on head side.
- Perineal Cut at foot side on the table & pair of lithotomy poles with straps
- Sliding footstep located at bottom center.
- Finish: Pre-treated & epoxy powder coated
- Overall Size (Approx.): 1870L x 560W x 820H mm.



EXAMINATION COUCH (PLAIN)

SMS-403

Description

- Fixed 50mm Foam Mattress with cover top in single section.
- Framework made of CRCA tubes mounted on PU stumps/shoes.
- One cabinet with one drawer.
- Finish: Pre-treated & epoxy powder coated
- Overall Approx. Size: 1830L x 610W x 810H mm



EXAMINATION TABLE

EXAMINATION TABLE MOTORISED (2 Section)

SMS-404

Description

- Fixed 50mm Foam Mattress with cover top in two section with single drawer.
- Framework made of CRCA tubes mounted on PU stumps/shoes.
- Back rest and Height adjustment by hand remote.
- Finish: Pre-treated & epoxy powder coated
- Overall Approx. Size: 1830L x 610W x 810H mm



GYNAE EXAMINATION TABLE SS (COUCH) Two Section

SMS-405

Description

- Fixed 50mm Foam Mattress with cover top in two section
- Framework made of Stainless Steel tubes mounted on PU stumps/shoes.
- Back rest adjustable by ratchet & Perineal Cut at footside with pair of lithotomy poles with straps
- Overall Approx. Size: 1830L x 560W x 810H mm

EXAMINATION TABLE SS (COUCH) Two Section

SMS-406

Description

- Fixed 50mm Foam Mattress with cover top in two section
- Framework made of Stainless Steel tubes mounted on PU stumps/shoes.
- Back rest adjustable by ratchet
- Overall Approx. Size: 1830L x 560W x 810H mm



EXAMINATION TABLE



EXAMINATION TABLE P/C (Two Section)

SMS-407

Description

- Fixed 50mm Foam Mattress with cover top in two section
- Framework made of CRCA tubes mounted on PU stumps/shoes.
- Back rest adjustable by ratchet.
- Finish: Pre-treated & epoxy powder coated
- Overall Approx. Size: 1830L x560W x 810H mm

GYNAE EXAMINATION TABLE (2 section)

SMS-408

Description

- Framework & 2 section top made of CRCA tubes & sheets mounted on PU stumps/shoes.
- Backrest adjustable by ratchet mechanism & Perineal Cut at footside with pair of lithotomy poles with straps
- Finish: Pre-treated & epoxy powder coated
- Overall Approx. Size: 1830L x560W x 810H mm

Optional

- Also available in stainless steel
- 2 section 50mm mattress with Raxin cover.



EXAMINATION TABLE (2 section)

SMS-409

Description

- Framework & 2 section top made of CRCA tubes & sheets mounted on PU stumps/shoes.
- Backrest adjustable by ratchet mechanism.
- Finish: Pre-treated & epoxy powder coated
- Overall Approx. Size: 1830L x560W x 810H mm

Optional

- Also available in stainless steel
- 2 section 50mm mattress with Raxin cover.



DELIVERY BED

LDR BED (Electric)

SMS-501



Description

- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Two individual Side Rails With Inbuilt Bed Angle Indicator
- Tubular framework mounted on 125mm (5") Central locking castors.
- Three-section top made of Bakelite sheet covered with 3 section Removable suitable Mattress.
- Hi – Low, Trendelenburg & Reverse Trendelenburg operated by Hand Remote .
- Leg section can smoothly slide under the main section.
- 1 Pair of lithotomy foam padded knee crutches.
- S.S Telescopic Heavy duty Drip Rod.

- Detachable S.S. Tray
- Provided with adjustable S.S. Grip Handles on both sides.
- Finish: Pretreated and epoxy powder coated

Dimensions (Approx.)

- Overall Size: 2000L x 950W mm
- Platform Size (full extended): 1900L x 660W
- Platform Size (Leg section Slide under): 1360L x 660W
- Minimum Height (without Mattress): 680mm
- Maximum Height (without Mattress): 960mm
- Backrest angle: 60°
- TB / RTB: 9° / 12°



DELIVERY BED

LDR BED (Hydraulic)

SMS-502

Description

- Polymer Moulded Head and Foot Panels easily removable without using any tools with corner buffer
- Polymer Moulded Two individual Side Rails With Inbuilt Bed Angle Indicator
- Tubular framework mounted on 125mm (5") Central locking castors.
- Three-section top made of Bakelite sheet covered with 3 section Removable suitable Mattress.
- Hi-Low adjustable by Hydraulic System (Foot paddle)
- Trendelenburg & Reverse Trendelenburg operated by Smooth crank system.
- Leg section can Smoothly slide under the main section.
- 1 Pair of lithotomy foam padded knee crutches.
- S.S Telescopic Heavy duty Drip Rod.
- Detachable S.S. waste collection Tray
- Provided with adjustable S.S. Grip Handles on both sides.
- Finish: Pretreated and epoxy powder coated.

Dimensions (Approx.)

- Overall Size: 2000L x 950W mm
- Platform Size (full extended): 1900L x 660W
- Platform Size (Leg section Slide under): 360L x 660W
- Minimum Height (without Mattress): 670mm
- Maximum Height (without Mattress): 920mm
- Backrest angle: 50°



OBSTETRIC LABOUR TABLE (Adjustable Height)

SMS-503

Description

- Stainless steel Head Side Rail easily drop down without using any tools with
- Stainless Steel Two drop down individual Side Rails
- Tubular framework mounted on 125mm (5") Central locking castors.
- Three-section top made of Stainless Steel sheet covered with 3 section Removable Mattress.
- Hi-Low adjustable by Hydraulic System (Foot paddle)
- Trendelenburg & Reverse Trendelenburg operated by Smooth crank system.
- Leg section can Smoothly slide under the main section.
- 1 Pair of lithotomy foam padded knee crutches.
- S.S Telescopic Heavy duty Drip Rod.
- Detachable S.S. waste collection Tray
- Provided with adjustable S.S. Grip Handles on both sides.
- Finish: Stainless Steel Matte & Glossy.

Dimensions

- Overall Size: 1900L x 900W mm
- Platform Size (full extended): 1800L x 660W mm
- Platform Size (Leg section Slide under): 1260L x 660W mm
- Adjustable Height (without Mattress): 670-920 mm
- Backrest angle: 50°
- TB / RTB: 9° / 12°

Note : Also Available in Epoxy Powder Coated (M.S)

DELIVERY BED

OBSTETRIC LABOUR TABLE (Fixed Height)

SMS-504



Description

- Stainless steel Head Side Rail easily drop down without using any tools with
- Stainless Steel Two drop down individual Side Rails
- Tubular framework mounted on PU Heavy Stump/ Shoes.
- Three-section top made of Stainless Steel sheet covered with 3 section Removable Mattress.

- Fix Height for Better Stability.
- Trendelenburg & Reverse Trendelenburg operated by Smooth crank system.
- Leg section can Smoothly slide under the main section.
- 1 Pair of lithotomy foam padded knee crutches.
- S.S Telescopic Heavy duty Drip Rod.
- Detachable S.S. waste collection Tray
- Provided with adjustable S.S. Grip Handles on both sides.
- Finish: Stainless Steel Matte & Glossy.

Dimensions

- Overall Size: 1900L x 900W x 720mm
- Platform Size (full extended): 1800L x 660W mm
- Platform Size (Leg section Slide under): 1260L x 660W mm
- Backrest angle: 50°
- TB / RTB: 9° / 12°

Note : – Also Available in Epoxy Powder Coated (M.S)



DELIVERY BED

OBSTETRIC DELIVERY TABLE (3 section)

SMS-505

Description

- Rectangular & square tube frame work mounted on PVC stumps.
- Three sections fixed cushioned top.
- Backrest, Seat & Leg section adjustable by ratchet mechanism.
- Provided with perineal 'U' cut.
- Provided with SS Waste collection bowls
- Single pair of Lithotomy crutches with Belt.
- The foot section of the bed can be adjustable- 0°-90°
- Knockdown construction.
- Finish: Pre treated & Epoxy Powder Coated.
- Overall Size (Approx.): 1830L x 610W x 810H mm



OBSTETRIC DELIVERY TABLE SS (3 section)

SMS-506

Description

- Rectangular & square tube frame work mounted on PVC stumps.
- Three sections fixed cushioned top.
- Backrest, Seat & Leg section adjustable by ratchet mechanism.
- Provided with perineal 'U' cut.
- Provided with SS Waste collection bowls
- Single pair of Lithotomy crutches with Belt.
- The foot section of the bed can be adjustable- 0°-90°
- Knockdown construction.
- Finish: all stainless steel matt/glossy
- Overall Size (Approx.): 1830L x 610W x 810H mm



WARD / OT EQUIPMENTS

CRASH CART S.S.

SMS-604

Description

- Framework made-up of Stainless Steel Tubes & Sheet.
- Provided with polystyrene boxes each with Lockable three drawers
- Provision of oxygen cylinder made of M.S.

- Provided with Sliding Massage board & which can be used as CPR board for cardiac patient
- Six colours removal Bins.
- 3 Stainless Steel Shelves for Different Uses.
- Provided with I.V. Rod & Corner Buffers.
- Cart mounted on 125 mm Diagonal locking castors.
- Overall Size (Approx.): 960L x 500 W x 1545H mm



M.S. Frame



S.S. Frame

MEDICINE TROLLEY (with 4 drawers)

SMS-605

Description

- Framework made-up of CRCA tube mounted on 10 cm (4") Diagonal Locking Castors.
- Two S.S. shelves with three side railing on top & Lower shelf
- Two M.S Drawer under the each shelf
- Finish: Pretreated & Epoxy Powder Coated.
- Overall Size (Approx.): 760L x 460W x 950H mm.



WARD / OT EQUIPMENTS

CRASH CART FULLY SS

SMS-606

Description

- Framework made-up of SS Tubes.
- SS storage cabinet with three drawers.
- Provision of M.S. oxygen cylinder holder
- Six colour removable bins.
- Provided with three SS shelves.
- Provided with SS IV rod and corner buffer and cardiac massage board.
- The cart is mounted on 125mm dia castors 2 with brakes.

Note : Also available in epoxy powder coated.



MULTI PURPOSE CRASH CART

SMS-607

Description

- Made of CRCA pipes/sheet
- Provided with six drawer of two different sizes
- First Three Drawers 520 x 410 x 50 mm
- Fourth & Fifth Drawer : 525 x 415 x 120mm
- Sixth Drawer : 535 x 425 x 185 mm
- Mounted on 125mm twin castors two with breaks.
- Accessories Oxygen cylinder Six coloured PVC Bins, Cardiac Massage Board, 2 Dustbin
- The complete metal components are pre treated and epoxy powder boated
- Overall Dimension : 1000L x 585W x 1125H mm



WARD / OT EQUIPMENTS

CASE HISTORY TROLLEY DOUBLE

SMS-608

Description:

- Framework & shelves made of ABS
- Mounted on 75mm castors
- 25 Case sheet Holders provision
- Approx. Sizes 340x400x850mm



CASE HISTORY TROLLEY SINGLE

SMS-609

Description:

- Framework & shelves made of ABS
- Mounted on 75mm castors
- 50 Case sheet Holders provision
- Approx. Sizes 680x400x850mm



WARD / OT EQUIPMENTS

INSTRUMENT TROLLEY (S.S) with Single Drawer

SMS-610

Description

- Framework made-up of S.S. tubes
- Two S.S. Shelves with S.S. railing on 3 sides
- The trolley mounted on 100 mm (4") diagonal locking Castors
- Finish: All stainless steel
- 610 L x 460 W x 815 H mm
- 760 L x 460 W x 815 H mm
- 915 L x 610 W x 815 H mm



INSTRUMENT TROLLEY (S.S) with Double Drawer

SMS-611

Description

- Framework made-up of S.S. tubes
- Two S.S. Shelves with S.S. railing on 3 sides
- The trolley mounted on 100 mm (4") diagonal locking Castors
- Finish: All stainless steel
- 610 L x 460 W x 815 H mm
- 760 L x 460 W x 815 H mm
- 915 L x 610 W x 815 H mm



CARE TROLLEY 3 SHELVES (S.S)

SMS-612

Description

- Framework made-up of S.S. tubes
- Three S.S. Shelves with central sliding drawer with 3 side railing.
- 6 rings on upper shelves
- The trolley mounted on 100 mm diagonal locking Castors
- Finish: All stainless steel
- 760 L x 460 W x 815 H mm
- 915 L x 610 W x 815 H mm

Optional : SS Dust Bin



WARD / OT EQUIPMENTS

LAPROSCOPIC / MONITOR TROLLEY

SMS-613

Description:

- Trolley made-up of CRCA sheet & Shelves & Single Drawer.
- Trolley has two handle for each of relocation.
- The trolley mounted on 125mm (5") Diagonal Locking Castors.
- Provided with Switch & Socket for medical equipment along with 1 Main MCB for protection
- Supplied with S.S. IV Rod and MS Oxygen Cylinder Cage.
- Finish: Pre-treated and Epoxy Powder Coated
- Overall Approx. Size: 700L x 620W x 1460H mm.



UTILITY TROLLEY

SMS-614

Description:

- Framework made of Tubular CRCA. pipe.
- Two Shelves made of Stainless steel.
- One Drawer made of CRCA.
- Trolley mounted on 100mm (4") Diagonal locking Castors
- Provided with S.S. I.V. Rod.
- Overall Size (Approx.): 520L x 450W x 940H mm.
- Top Size (Approx.): 450L x 450W mm.



ECG TROLLEY

SMS-615

Description:

- Framework made of Tubular CRCA. pipe.
- Two Shelves made of Stainless steel.
- One Drawer made of CRCA.
- Trolley mounted on 100mm (4") Diagonal locking Castors
- Provided with S.S. I.V. Rod.
- Overall Size (Approx.): 520L x 450W x 940H mm.
- Top Size (Approx.): 450L x 450W mm.



WARD / OT EQUIPMENTS

DRESSING TROLLEY

SMS-616

Description

- Framework made-up of S.S. tubes
- Two S.S. Shelves with S.S. railing on 3 sides
- The trolley fitted S.S. bucket and S.S bowl
- The trolley mounted on 100 mm diagonal locking Castors
- Bowl & Bucket Holding cage on One Side in Knockdown position.
- Finish: All stainless steel
- Overall Size (Approx.): 1010L x 460W x 810H mm.



INSTRUMENT TROLLEY (S.S)

SMS-617

Description

- Framework made-up of S.S. tubes
- Two S.S. Shelves with S.S. railing on 3 sides
- The trolley mounted on 100 mm diagonal locking Castors
- Finish: All stainless steel
- 610 L x 460 W x 815 H mm
- 760 L x 460 W x 815 H mm
- 915 L x 610 W x 815 H mm

INSTRUMENT TROLLEY 3 SHELVES (S.S)

SMS-618

Description

- Framework made-up of S.S. tubes
- Three S.S. Shelves with S.S. railing on 3 sides
- The trolley mounted on 100 mm diagonal locking Castors
- Finish: All stainless steel
- 760 L x 460 W x 815 H mm
- 915 L x 460 W x 815 H mm
- 915 L x 610 W x 815 H mm



WARD / OT EQUIPMENTS

MAYO'S TROLLEY (Double Gear)

SMS-619

Description:

- Both side Frame work made of S.S. rectuangler / Square tubes
- Height adjustable by both side gear handle
- Tray made of 0.9 mm thick S.S. 304 Grade sheet
- Table mounted on 100 mm Diagonal locking castors
- Finish: S.S. Frame with S.S. Removable tray
- Overall Size (Approx.): 1220L x 520W x 910H mm
- Tray Size (Approx.): 925L x 545W mm



MAYO'S TROLLEY (Gear Single)

SMS-620

Description:

- Frame work made of S.S. rectuangler / Square tubes
- Height adjustable by Single gear handle
- Table mounted on 50 mm Diagonal locking castors
- Finish: S.S. Frame with S.S. Removable tray
- Overall Size (Approx.): 600L x 400W x 900 – 1200H mm.
- Tray Size (Approx.): 625L x 440W mm

Note: Also available M.S Framework S.S Tray.

MAYO'S TROLLEY (Knob)

SMS-621

Description:

- Frame work made of S.S. rectuangler / Square tubes
- Height adjustable by manual Knob.
- Table mounted on 50 mm Diagonal locking castors
- Finish: S.S. Frame with S.S. Removable tray
- Overall Size (Approx.): 600L x 400W x 900 – 1200H mm.
- Tray Size (Approx.): 625L x 440W mm

Note: Also available M.S Framework S.S Tray.



WARD / OT EQUIPMENTS



SURGEON CHAIR

SMS-701

Description

- Cushioned top with 5 leg base mounted on 5cms castors.
- Adjustable height with pneumatic pump
- Provided with back support and tubular ring for foot support



SURGEON STOOL

SMS-702

Description

- Cushioned top with 5 leg base mounted on 5cms castors.
- Adjustable height with pneumatic pump
- Tubular ring for foot support



REVOLVING STOOL (All SS)

SMS-703

Description

- Four SS. Legs base fitted with PVC Stumps/Shoes.
- Height adjustable by mechanical screw.
- SS Revolving top.
- Finish: Stainless Steel matte & glossy
- Overall Size (Approx.): 480L x 480W x 500H mm.

Also available in P.C Frame work with SS Top

REVOLVING STOOL (Cushioned Top)

SMS-704

Description

- Four legs S.S. base fitted with PU Stumps/Shoes
- Height adjustable by mechanical screw
- Moulded Cushioned Revolving top.
- Finish: SS glossy finish.
- Overall Size (Approx.): 510Lx 510W x 540H mm.

Also available in P.C Frame work

WARD / OT EQUIPMENTS



VISITOR STOOL

SMS-705

Description

- Frame & Top made of SS. Sheet & SS Tube
- Legs fitted with PU Stumps/Shoes.
- Finish: Stainless Steel matte & glossy finish.
- Overall Size (Approx.): 350L x 350W x 490H mm

Note: Also available CRCA Framework.



TROLLEY FOR OXYGEN CYLINDER

SMS-706

Description

- Framework made of CRCA tube mounted on 100 mm castors at rear.
- Finish: Pre-treated and Epoxy powder coated
- For D type cylinder (Big Cylinder)
- For B type cylinder (Small Cylinder)



WASH BASIN STAND (Double)

SMS-707

Description:

- Basin holder Ring & Vertical Pipe made of CRCA tube
- Mounted on MS Base with 50mm dia castor.
- 2 Stainless Steel basin dia 16"
- Finish: Epoxy Powder Coated.



WASH BASIN STAND (Single)

SMS-708

Description:

- Basin holder Ring & Vertical Pipe made of CRCA tube
- Mounted on MS Base with 50mm dia castor.
- Single Stainless Steel basin dia 16"
- Finish: Epoxy Powder Coated

WARD / OT EQUIPMENTS



FOOT STEP SINGLE P/C

SMS-709

Description

- Made from S.S. Sheet mounted on PVC stumps
- Dimensions : 18"L x 9"W x 9"H



FOOT STEP SINGLE SS

SMS-710

Description

- Made from S.S. Sheet mounted on SS stumps
- Dimensions : 18"L x 9"W x 9"H



FOOT STEP SINGLE MOLDED RUBBER

SMS-711

Description

- Stainless steel tubular legs fitted with PVC stump
- Non-skid molded rubber matting steps top
- Knock down construction
- Dimensions : 400W x 280D x 225H mm

Also available in Epoxy Powder coating



FOOT STEP DOUBLE SS

SMS-712

Description

- Made from S.S. Sheet mounted on PVC stumps
- Dimensions : 18"L x 18"W x 18"H
- First Step Height : 230mm
- Second Step Height : 450mm

Also available in CRCA sheet Epoxy Powder coated



FOOT STEP DOUBLE SS MOLDED RUBBER

SMS-713

Description

- Stainless steel tubular legs fitted with PVC stump
- Non-skid molded rubber matting steps top
- Knock down construction
- Dimensions : 525W x 400D x 350H mm



FOOT STEP DOUBLE P/C MOLDED RUBBER

SMS-714

Description

- CRCA tubular legs fitted with PVC stump
- Non-skid molded rubber matting steps top
- Pre treated and Epoxy Powder Coated
- Knock down construction
- Dimensions : 525W x 400D x 350H mm

WARD / OT EQUIPMENTS

BED SIDE SCREEN SS (3 panels)

SMS-715

Description

- Framework made of S.S. Tubular in three section.
- Legs mounted on 50mm dia castors with hooks and spring wire.
- Plastic washable curtain in blue colour.
- Finish: Stainless Steel Glossy
- Overall Size: 1727H x 1829L mm



MS/P.C.

BED SIDE SCREEN SS (4 panels)

SMS-716

Description

- Framework made of S.S. Tubular in Four section.
- Legs mounted on 50mm dia castors with hooks and spring wire.
- Plastic washable curtain in blue colour.
- Finish: Stainless Steel Glossy
- Overall Size: 1727H x 2438L mm



MS/P.C.



WARD / OT EQUIPMENTS



DRESSING DRUMS (SEAMLESS SS)

SMS-717

Available Sizes:

6"x4", 6"x6", 9"x5", 11"x5", 11"x9",
14"x9", 15"x12"



INSTRUMENT TRAY (SEAMLESS SS)

SMS-718

Available Sizes:

8"x6", 9"x6", 10"x8", 11"x7",
12"x8", 12"x10", 15"x12", 18"x12"



CATHETER TRAY

SMS-719

Available Sizes:

8"x3", 16"x4"



BED PAN (SS)

SMS-720

Available for Male / Female



URINAL POT (SS)

SMS-721

Available for Male / Female



BASIN/BOWL (SS)

SMS-722

Available Sizes: 12", 16"



LOTION BOWL (SS)

SMS-723



KIDNEY TRAY (SEAMLESS SS)

SMS-724

Available Sizes:

6", 8", 10", 12"



ENEMA CANE (SS)

SMS-725

WARD / OT EQUIPMENTS



KICK BUCKET

SMS-726

Description

- Tubler framework mounted on 50cms castors.
- Stainless Steel Bucket Cap. 12Ltrs.
- Pre Treated and Epoxy powder coated.
- Size : H12" x 14" Dia
- Also available in Stainless steel.



KICK BUCKET

SMS-727

Description

- SS Bucket with stainless steel frame with castors.
- Dimension : 35 D x 30H cm



DUST BIN

SMS-728

Description

- Foot pedal Operated.
- Stainless Steel construction.
- Removable plastic bin
- Base fitted with Anti-Skid Plastic.
- Capacity : 10/15/20/30 Ltrs

BIO MEDICAL WASTE BIN (Only Bin)

SMS-729

Description

- Bins made up of plastic.
- Bins are foot operated
- Available in red / blue / black/ yellow/ white / green colours.
- Available sizes: 20 ltrs, 30 ltrs, 80 ltrs with 2 pcs rear wheels, 100-120 ltrs with 2 pcs rear wheels but without foot operated.



BIO MEDICAL WASTE BIN (SETS)

SMS-730

Description

- Bins made up of plastic.
- Framework made up of SS tubes.
- Bins are foot operated
- Sets of 4 bins / 3 bins / 2 bins can be altered as required, standard is 4 bins in 1 set
- Available in red / blue / black/ yellow/ white / green colours



WARD / OT EQUIPMENTS

SOILED LINEN TROLLEY WITH PLASTIC BUCKET

SMS-736

Description

- Framework & Push handle made of SS tubular
- Mount on 100mm Diagonal castors
- 100ltrs. Capacity Premium quality washable plastic bucket.
- Finish; Stainless steel glossy .
- Aprox. Sizes 650x530x930mm.

Note: also available Epoxy Powder Coated framework



SOILED LINEN TROLLEY WITH CANVAS BAG

SMS-737

Description

- Framework made of SS tubular
- Mount on 100mm Diagonal castors
- Superior quality of canvas bag.
- Finish; Stainless steel glossy.
- Aprox. Sizes 520x515x910mm

SALINE STANDS

SMS-738

Description

- SS pipe & rod (4 hooks) with ss base
- SS pipe & rod (2 hooks) with ss base
- SS pipe & rod (4 hooks) with plastic base
- SS pipe & rod (2 hooks) with plastic base
- MS pipe & SS rod with MS base
- MS pipe & SS rod with plastic base



ACCESSORIES

MATTRESS FOR HOSPITAL BEDS

SMS-731

Description


- 4 Section PU Foam 4" Mattress with Premium Zipper Rexine Cover for ICU & Fowler Bed



SMS-732

Description


- 2 Section PU Foam 4" Mattress with Premium Zipper Rexine Cover for Semi Fowler Bed



SMS-733

Description

- Single Section PU Foam 4" Mattress with Premium Zipper Rexine Cover for Plain & Pediatric Bed



BED SIDE RAILINGS

SMS-734



Swing Down Railing (SS & MS)

SMS-735



Imported Collapsible Railing

SMS-739



ABS Railing

SMS-740



Collapsible Railing (SS & MS)

DIALYSIS BED / CHAIR

DIALYSIS BED (MOTORIED)

SMS-801

Description

- Framework made of CRCA Tubular & mounted on 125mm Diagonal Locking Heavy duty castors.
- Three-section Platform top with Adjustable height with high density Moulded PU Mattress.
- Operated: Backrest, Kneerest & Hi – Low by hand remote
- Single Button Chair Position.
- Easy detachable head & foot panel (foam padded panel)
- Mechanism covered by heavy duty CRCA base.
- Telescopic heavy I.V. Rod with 4 hooks
- Foam padded adjustable handrest on both sides
- Finish: Pretreated and epoxy powder coated.

Dimensions (Approx.):

- Overall Size: 1900L x 750W x 670-970 mm.
- Backrest : 85°
- Kneerest : 80° (Down)



BLOOD TRANSFUSION CHAIR

SMS-802

Description

- Framework made of CRCA tube.
- The fabric is pressed PU artificial leather
- Three angles adjustable
- Single side basket & IV Rod .
- Finish: Epoxy Powder coated
- Open Size: 700 x 770 x 1080 mm



BLOOD SAMPLER CHAIR

SMS-803

Description

- Framework made of CRCA tube.
- The back and seat padded with PU Foam covered from good quality leatherised Rexin cover.
- Finish: Epoxy Powder coated
- Size: Seat length 520L x 450W x Backrest Height 850 mm



GENERAL EQUIPMENTS

WHEELCHAIR NON-FOLDING

SMS-804

Features and Benefits

- Non folding 17" seat and fixed cushioned back rest, fixed arm rest, 5" front castors, 26" rear wheel with brakes.
- Pre-treated Epoxy powder coated
- Dimension : 39"L x 25"W x 36"H

Note : Also available in Stainless Steel



WHEELCHAIR FOLDING

SMS-805

Features and Benefits

- Foldable wheelchair 17" square seat, fixed arm-rest, 26" rear wheel, 6" front castors.
- Pre-treated powder coated
- Dimension : 42"L x 26"W x 36"H

ONE USER SCRUB STATION ELBOW OPERATED

SMS-806

Features and Benefits

- Manufactured with heavy gauge stain finish sheets for durability.
- Size : 900L x 540W x 840H mm



TWO USERS SCRUB STATION FOOT OPERATED

SMS-807

Features and Benefits

- Manufactured with heavy gauge stain finish sheets for durability.
- Size : 1450L x 540W x 840H mm

THREE USERS SCRUB STATION SENSOR CONTROL

SMS-808

Features and Benefits

- Manufactured with heavy gauge stain finish sheets for durability.
- Size : 1950L x 540W x 840H mm



OUR PLANTS





Gr. Floor, Plot No. 326, Green Park, Sr. No. 96/1/1, Near HP Petrol Pump,
Lohegaon, Pune, Maharashtra - 411047

Mob. : +91 9140230733, 7070708204

E-mail : shreyamedicalsystem2010@gmail.com

Website : www.shreyamedicalsystem.com, www.shreyamodularot.co.in